

# Grid Pad Specifications January 2017



Go 8



Go 10



Pro 11



Pro 13



Pro 18



Eye 11



Eye 13



Eye 18

Operating system	Windows 10	Windows 10	Windows 10	Windows 10	Windows 10	Windows 10	Windows 10	Windows 10
Software	Grid 3	Grid 3	Grid 3	Grid 3	Grid 3	Grid 3	Grid 3	Grid 3
Processor	Intel® Atom™ processor Z3795 1.6GHz to 2.4 GHz Burst	Intel® Atom™ processor Z3795 1.6GHz to 2.4 GHz Burst	Intel® Core™ processor m5-6Y54 1 GHz - 2.7 GHz Burst	Intel Celeron 2957U 1.40GHz	Intel® Pentium® processor J3710 1.6 GHz - 2.6 GHz Burst	Intel® Core™ processor m5-6Y54 1 GHz - 2.7 GHz Burst	Intel Celeron 2957U 1.40GHz	Intel® Pentium® processor J3710 1.6 GHz - 2.6 GHz Burst
Memory	4GB	4GB	4GB	8GB	8GB	4GB	8GB	8GB
Hard drive	64GB	64GB	128GB	256GB SSD	256GB SSD	128GB	256GB SSD	256GB SSD
Video card	Intel® HD Graphics	Intel® HD Graphics	Intel® HD Graphics	Intel® HD Graphics	Intel® HD Graphics	Intel® HD Graphics	Intel® HD Graphics	Intel® HD Graphics
Dimensions	227.8 x 12.5 x 150.1mm 8.97x 0.49 x 5.91 inch	269 x 180 x 12.5 mm 10.6 x 7.09 x 0.49 inch	310 x 200 x 50 mm 12.2 x 7.8 x 2 inch	344 x 216 x 57 mm 13.5 x 8.5 x 2.2 inch	440 x 325 x 50 mm 17.3 x 12.8 x 2 inch	310 x 200 x 50 mm 12.2 x 7.8 x 2 inch	344 x 216 x 57 mm 13.5 x 8.5 x 2.2 inch	440 x 325 x 50 mm 17.3 x 12.8 x 2 inch
Weight	0.54 kg 1.19 lbs	0.75 kg 1.65 lbs	1 kg 2.2 lbs	2.68 kg 5.9 lbs	2.75 kg 6 lbs	1.7 kg* 3.7 lbs*	3.1 kg* 6.8lbs*	3.25 kg* 7.17 lbs*
Battery life	Up to 5hrs	Up to 6hrs	Up to 6hrs	Up to 8hrs, hot swap batteries provided	Up to 4hrs 40 mins	Up to 2 hrs 30 mins with eye gaze	Up to 5 hrs with eye gaze, hot swap available	Up to 2 hrs with eye gaze
Tablet camera	2MP front, 8MP rear	2MP front, 8MP rear	2MP front, 5MP rear	2MP front, 5MP rear	2MP front	2MP front, 5MP rear	2MP front, 5MP rear	2MP front
Screen	8.3" 1920 x 1200 HD display	10.1" 1920 x 1200 HD display	11.6" 1080 HD display	13.3" LED 1080 Full HD display	17.3" 1080 full HD display	11.6" 1080 HD display	13.3" LED 1080 Full HD Display	17.3" 1080 full HD display
Eye gaze camera	n/a	n/a	n/a	n/a	n/a	Alea IntelliGaze, MyGaze	Alea IntelliGaze, MyGaze	Alea IntelliGaze, MyGaze
Access	Touchscreen, USB input	Touchscreen, USB input	Touchscreen, switch, USB input	Touchscreen, switch, USB input	Touchscreen, switch, USB input	Touchscreen, switch, USB input	Touchscreen, switch, USB input	Touchscreen, switch, USB input
Connectivity	Wi-fi, bluetooth	Wi-fi, bluetooth	Wi-fi, bluetooth	Wi-fi, bluetooth	Wi-fi, bluetooth	Wi-fi, bluetooth	Wi-fi, bluetooth	Wi-fi, bluetooth
Ports	1 x USB2.0, speaker out, micro SD card slot (Micro USB charging only)	1 x USB2.0, speaker out, micro SD card slot (Micro USB charging only)	3 x USB2.0, speaker out, micro SD, micro HDMI, auditory feedback port	2 x USB 2.0, Auditory feedback, speaker out, microphone	2 x USB 2.0, 1 x USB3.0, Speaker Out, Auditory feedback, headphone, SD card reader, auditory feedback socket, HDMI	3 x USB2.0, speaker out, micro SD, micro HDMI, auditory feedback port	2 x USB 2.0, speaker out, auditory feedback, reader, speaker out, microphone	2 x USB 2.0, 1 x USB3.0, speaker out, micro SD card reader, auditory feedback socket, HDMI
Buttons	Power, volume	Power, volume	Power, volume	Power	Power, home, volume	Power, volume	Power	Power, home, volume
Switch access	n/a	n/a	2 x 3.5mm jacks, radio switch (additional extra)	2 x 3.5mm jacks, radio switch (additional extra)	2 x 3.5mm jacks, radio switch (additional extra)	2 x 3.5mm jacks, radio switch (additional extra)	2 x 3.5mm jacks, radio switch (additional extra)	2 x 3.5mm jacks, radio switch (additional extra)
Amplified speaker	External bluetooth speaker	External bluetooth speaker	2 x built-in high quality speakers	2 x built-in high quality speakers	2 x built-in high quality speakers	2 x built-in high quality speakers	2 x built-in high quality speakers	2 x built-in high quality speakers
Case	Integrated rugged protection with hand and shoulder strap (optional Tough Case)	Integrated rugged protection with hand and shoulder strap (optional Tough Case)	Integrated Tough Case	n/a	n/a	Integrated Tough Case	n/a	n/a
Mounting options	n/a	n/a	Rehadapt, Daessy and built in Tablet Stand	Rehadapt, Daessy and built in Tablet Stand	Rehadapt and Daessy	Rehadapt, Daessy and built in Tablet Stand	Rehadapt, Daessy and built in Tablet Stand	Rehadapt and Daessy
Environment control	n/a	n/a	Servus inside - Gewa, Eldat or Z-wave, radio switch receiver	Servus inside - Gewa, Eldat or Z-wave, radio switch receiver	Servus inside - Gewa, Eldat or Z-wave, radio switch receiver	Servus inside - Gewa, Eldat or Z-wave, radio switch receiver	Servus inside - Gewa, Eldat or Z-wave, radio switch receiver	Servus inside - Gewa, Eldat or Z-wave, radio switch receiver
Warranty	2 year free warranty	2 year free warranty	2 year free warranty	2 year free warranty	2 year free warranty	2 year free warranty	2 year free warranty	2 year free warranty

\*System with myGaze Assistive 2 eye gaze camera