



Device	Windows based AAC devices, laptop or PC	iOS 11 and above including iPad mini, Air and Pro
Operating system	Windows 7, 8, or 10	iOS 12 and above
Versions	<p>Grid 3 60 day free trial</p> <p>Grid 3 Education Licence Ideal for schools and organisations with managed networks</p>	<p>Grid for iPad - AAC 30 days free trial, monthly subscription or complete purchase</p> <p>Grid for iPad - VPP Complete purchase only, ideal for schools and organisations with MDM</p>
Access	All Access	Touch and Switch
Key features	<ul style="list-style-type: none"> Chat history SwiftType Prediction Message banking Remote editing Online Grids Accessible apps Computer control Environment control 	<ul style="list-style-type: none"> Chat history SwiftType Prediction Message banking Remote editing Online Grids
Symbol communication content	<ul style="list-style-type: none"> Super Core Voco Chat Word Power PODD (separate purchase) 	<ul style="list-style-type: none"> Super Core Voco Chat PODD (separate purchase)
Text communication content	<ul style="list-style-type: none"> Text Talker Fast Talker Alpha Core 	<ul style="list-style-type: none"> Text Talker Alpha Core
Interactive Learning content	<ul style="list-style-type: none"> Cause and Effect Choices Challenges Visual Scenes First Words 	<ul style="list-style-type: none"> Cause and Effect Choices Challenges Visual Scenes First Words
Symbol libraries	Widget, SymbolStix, PCS, PCS High Contrast, Bliss, Metacom, Snaps Photos, Tawasol, See Sense, Beta	Widget, SymbolStix, PCS, PCS High Contrast, Snaps Photos, Tawasol, Metacom
Voices	Over 100 voices in 30+ languages Support for Accapela, Nuance, SAPI, CereVoice Me, SpeakUnique	Over 100 voices in 30+ languages Support for Accapela, Nuance, iOS, SpeakUnique