



Power Pad is our most powerful communication aid yet.

This all-in-one AAC device is designed for people with complex access needs, and will work with any alternative access method.

Operating system	Windows 10 Pro
Software	Grid 3
Processor	Intel® Core™ i5-6200U Skylake processor, 3MB, 2.1GHz up to 2.7 GHz
Memory / RAM	8GB DDR3
Hard drive	SSD SATA III, 128GB
Dimensions	355 x 266 x 70 mm 14 x 10.5 x 2.7 inch
Weight	3.4 kg 7.5 lbs
Battery life	15 hours typical use 9 hours with eye gaze Additional Hot swappable batteries available to purchase - 6 hours each
Tablet camera	2MP front, 5MP rear
Display	13.3" Gorilla Glass Intel HD 1920 X 1080p
Eye gaze	Eyetech TM5 Mini Head motion box (31.5 X 22.5 X 35cm) Eyewear compatibility
Access	2 x 3.5mm switch ports, 3.5mm headphone / mic combo port, touch screen, Radio switch input, magnetic Smartnav mount, Keyguard compatible, 2 x USB 2.0 ports, Radio switch receiver (RF)
Environment control	GEWA IR, Easywave / Z-Wave Environment control
Connectivity	Dual band Wireless (11ac/abgn), Bluetooth 4.0 LE
Mounting plate	Daessy and Rehadapt
Other features	Variable angle stand
Warranty	3 years with SmartCare (UK and Ireland only)