



Grid Pad 12 is designed for communication all day, every day. A long 15-hour battery life, rugged design and multiple access options ensure you always have your voice.

Operating system	Windows 10 Pro
Software	Grid 3
Processor	Intel® Core™ i5-7Y54 Processor, up to 3.20 GHz
Memory / RAM	8GB DDR3
Hard drive	SSD SATA III, 256GB
Display	12.5" 1920 X 1080p full HD with hardened glass
Weight	2.7kg 5.9 lbs
Battery life	15 hours typical use 10 hours with eye gaze
Charge time	10 hours
Tablet camera	2MP front, 5MP rear
Dimensions	355 x 223 x 45 mm 14 x 8.8 x 1.8 inch
Eye gaze	Optional Grid Pad camera unit (Alea IntelliGaze, EyeTech TM5 Mini or Irisbond Duo) Eye gaze mounting plates also available (PCEye Mini, LC Edge and MyGaze)
Access	2 x 3.5mm switch ports, 3.5mm headphone / mic combo port, touch screen, magnetic Keyguard compatible, 1 x USB 2.0 port, 1 x USB3.0 port (Eye Gaze uses the USB 3.0 port)
Environment control	GEWA Infrared, EasyWave / Z-Wave radio
Connectivity	Dual band Wireless (AC 8265), Bluetooth 4.2 LE
Mounting plate	Daessy and Rehadapt
Warranty	3 years SmartCare (UK only)
Other features	Built in stand and Remote Power Button (batteries included)