



The Smartbox Essentials range brings together the hardware, software, accessories and resources to help you get started with AAC.

iPad mini

iPad

Hardware	Apple iPad mini 5th Generation 2019, 64GB	Apple iPad 7th Generation 2019, 32GB
Operating system	iOS 11	iOS 11
Software	Grid for iPad	Grid for iPad
Case	Soft foam case (orange) or Tough shell case	Soft foam case (blue) or Tough shell case
Performance	64-bit A9 chip, Solid state hard drive 64GB	64-bit A9 chip, Solid state hard drive 32GB
Connectivity	Dual band Wireless Bluetooth 4.2	Dual band Wireless Bluetooth 4.2
Display	7.9" Durable Retina display touchscreen	9.7" Durable Retina display touchscreen
Battery life	Up to 10 hours	Up to 10 hours
Amplification	Smartbox bluetooth speaker	Smartbox bluetooth speaker
Weight	Device with Soft foam case: 444g Device with Tough shell case: 527g	Device with Soft foam case: 730g Device with Tough shell case: 910g
Ports	1 X charging port , 3.5mm headphone jack / microphone port	1 X charging port , 3.5mm headphone jack / microphone port
Camera	Front and rear camera	Front and rear camera
Accessories	Carry case, extra charger, AAC classroom pack, Strap (Tough shell case only)	Carry case, extra charger, AAC classroom pack, Strap (Tough shell case only)