

Grid Pad comparison

| | Touch Pad | Touch Pad Upgrade | Grid Pad 8 | Grid Pad 10 | Grid Pad 12 |
|----------------------|---------------------------------|-------------------------------|----------------------------------|----------------------------------|------------------------------|
| HD display | 10.5" | 10.5" | 8" | 10" | 12" |
| Weight | 1062g | 1062g | 1.2kg | 1.4kg | 2.7kg |
| Battery life | 9.5 hours | 7 hours | 4 hours | 7 hours | 15 hours |
| Processor | Intel®Pentium Gold 4425Y 1.7GHz | Intel® Core™ m3-8100Y 1.1 GHz | Atom z3795 1.6GHz (2.4GHz Burst) | Atom z3795 1.6GHz (2.4GHz Burst) | Intel® Core™ i5-7Y54 3.2 GHz |
| Memory / RAM | 4GB | 8GB | 4GB | 4GB | 8GB |
| USB ports | 1 X USB C | 1 X USB C | 2 X USB 2.0 | 2 X USB 2.0 | 1 X USB 2.0, 1 X USB 3.0 |
| Windows 10 | ✓ | ✓ | ✓ | ✓ | ✓ |
| Grid 3 | ✓ | ✓ | ✓ | ✓ | ✓ |
| Tablet camera | ✓ | ✓ | ✓ | ✓ | ✓ |
| Remote power button | | | | | ✓ |
| Eye gaze | | | | | ✓ |
| Switch | | | ✓ | ✓ | ✓ |
| Headphone port | ✓ | ✓ | ✓ | ✓ | ✓ |
| Environment control | | | ✓ | ✓ | ✓ |
| Mountable | ✓ | ✓ | ✓ | ✓ | ✓ |
| Portable | ✓ | ✓ | ✓ | ✓ | |
| Strap and tough case | ✓ | ✓ | ✓ | ✓ | |