

Touch Pad



Touch Pad is a lightweight device designed for touch access. It's perfect for users who need a powerful and portable device with a long battery life and loud speaker.

	Touch Pad	Touch Pad Upgrade
Operating system	Windows 10 Pro	Windows 10 Pro
Software	Grid 3	Grid 3
Processor	Intel® Pentium® Gold Processor 4425Y 1.70GHz, 2 MB Intel® Smart Cache	Intel® Core™ m3-8100Y Processor 4M Cache, 1.1 GHz base, 3.40 GHz burst
Memory / RAM	4GB	8GB
Hard drive	64GB	128GB
Display	10.5" 1920 x 1280	10.5" 1920 x 1280
Access	Touch screen	Touch screen
Weight	1062g	1062g
Battery life	9.5 hours	7 hours
Charge time	Up to 3 hours	Up to 3 hours
Tablet camera	5MP front, 8MP rear	5MP front, 8MP rear
Dimensions	260 x 215 x 60 mm	260 x 215 x 60 mm
Connectivity	Wireless (IEEE 802.11 a/b/g/n), Bluetooth 5.0	Wireless (IEEE 802.11 a/b/g/n), Bluetooth 5.0
Mounting plate	Rehadapt and Daessy	Rehadapt and Daessy
Other features	Screen protector Handle Strap Speaker Tough case Keypad ready USB C to A adapter	Screen protector Handle Strap Speaker Tough case Keypad ready USB C to A adapter