

Grid Pad comparison

	Touch Pad	Grid Pad 8	Grid Pad 10	Grid Pad 12
HD display	10"	8"	10"	12"
Weight	1kg	1.2kg	1.4kg	2.7kg
Battery life	8 hours	4 hours	7 hours	15 hours
Processor	Intel®Pentium Gold 4415Y 1.6GHz	Atom z3795 1.6GHz (2.4GHz Burst)	Atom z3795 1.6GHz (2.4GHz Burst)	Intel® Core™ i5-7Y54 3.2 GHz
Memory / RAM	4GB	4GB	4GB	8GB
USB ports	1 X USB C	2 X USB 2.0	2 X USB 2.0	1 X USB 2.0, 1 X USB 3.0
Windows 10	✓	✓	✓	✓
Grid 3	✓	✓	✓	✓
Tablet camera	✓	✓	✓	✓
Remote power button				✓
Eye gaze				✓
Switch		✓	✓	✓
Headphone port	✓	✓	✓	✓
Environment control		✓	✓	✓
Mountable	✓	✓	✓	✓
Portable	✓	✓	✓	
Shoulder strap and tough case	✓	✓	✓	