

Touch Pad



Surface Touch is a lightweight device designed for touch access. It's perfect for users who need a powerful and portable device with a long battery life and loud speaker.

Operating system	Windows 10 Pro
Software	Grid 3
Processor	Intel®Pentium Gold 4415Y 1.6GHz
Memory / RAM	4GB
Hard drive	64GB
Display	10" 1800 x 1200
Access	Touch screen
Weight	1Kg
Battery life	8 hours
Tablet camera	5MP front, 8MP rear
Dimensions	260 x 215 x 60 mm
Connectivity	Wireless (IEEE 802.11 a/b/g/n/ac), Bluetooth 4.1
Mounting plate	Rehadapt and Daessy
Other features	Handle Strap Speaker Tough case Keypad ready USB C to A adapter