



Grid Pad 15 is designed for communication all day, every day. A long 14-hour battery life, rugged design and multiple access options ensure you always have your voice.

Operating system	Windows 10 Pro
Software	Grid 3
Processor	Intel® Core™ i5-7Y54 Processor (1.2GHz, 3.2GHz Burst)
Memory / RAM	8GB
Hard drive	SSD SATA III, 256GB
Display	15.6" 1920 x 1080 with toughened glass and anti-glare finish
Second screen	4.6" 800 x 300 with toughened glass and anti-glare finish
Weight	3kg 6.6 lbs
Battery life	14 hours typical use 9 hours with eye gaze
Charge time	10 hours
Tablet camera	2MP front, 5MP rear
Dimensions	408 x 259 x 42 mm 21 x 10 x 1.7"
Eye gaze	Optional Grid Pad camera unit (Alea IntelliGaze, EyeTech TM5) Eye gaze mounting plates also available (PCEye Mini, Eyegaze Edge and MyGaze)
Access	2 x 3.5mm switch ports, 3.5mm headphone / mic combo port, touch screen, magnetic keyguard compatible, 1 x USB 2.0 port, 1 x USB 3.0 port (Eye gaze uses the USB 3.0 port)
Environment control	GEWA Infrared, EasyWave / Z-Wave radio
Connectivity	Dual band wireless (AC 8265), Bluetooth 4.2 LE
Mounting plate	Daessy and Rehadapt
Warranty	3 years SmartCare (UK only)
Other features	Antimicrobial casing and screen, built in stand and Remote Power Button