



Grid Pad 10s

Windows 10 Pro

Intel® Core™ i5-7Y54 processor

8 GB RAM

256 GB SSD

10.1" 1920 x 1200 daylight bright display

3.9" 480 x 128 second screen

1.45kg / 3.2lb

7h battery life (typical use)

5h charge time

Dimensions: 260 x 189 x 60 mm

2 x 3.5mm switch ports

2 x USB 3.0 ports

3.5mm headphone / mic port

GEWA Infrared Environment Control

Dual band wireless (AC 8265), Bluetooth 4.2 LE

Daessy and Rehadapt mounting

Antimicrobial

Remote Power Button

IP54 (dust and splash resistant)



Grid Pad 12

Windows 10 Pro

Intel® Core™ i5-7Y54 processor

8 GB RAM

128 GB SSD

12.5" 1920 x 1080 display

Optional second screen accessory

2.7kg / 6lb

15h battery life (typical use)

10h charge time

Dimensions: 355 x 223 x 45 mm

2 x 3.5mm switch ports

1 x USB 2.0 port, 1 x USB 3.0 port

3.5mm headphone / mic port

GEWA Infrared and Radio Environment Control

Dual band wireless (AC 8265), Bluetooth 4.2 LE

Daessy and Rehadapt mounting

Antimicrobial

Remote Power Button

Optional eye gaze module



Grid Pad 15

Windows 10 Pro

Intel® Core™ i5-7Y54 processor

8 GB RAM

256 GB SSD

15.6" 1920 x 1080 display

4.6" 800 x 300 second screen

3kg / 6.6lb

14h battery life (typical use)

10h charge time

Dimensions: 408 x 259 x 42 mm

2 x 3.5mm switch ports

1 x USB 2.0 port, 1 x USB 3.0 port

3.5mm headphone / mic port

GEWA Infrared and Radio Environment Control

Dual band wireless (AC 8265), Bluetooth 4.2 LE

Daessy and Rehadapt mounting

Antimicrobial

Remote Power Button

Optional eye gaze module

