

Touch Pad



Touch Pad is a lightweight device designed for touch access. It's perfect for users who need a powerful and portable device with a long battery life and loud speaker.

	Touch Pad	Touch Pad Upgrade
Operating system	Windows 10 Pro	Windows 10 Pro
Software	Grid 3	Grid 3
Processor	Intel® Pentium® Gold 6500Y 3.40GHz	Intel® Core™ i3-10100Y, 3.9GHz
Memory / RAM	4GB	8GB
Hard drive	64GB	128GB
Display	10.5" 1920 x 1280	10.5" 1920 x 1280
Access	Touch screen	Touch screen
Weight	1062g	1062g
Battery life	8.5 hours	7 hours
Charge time	Up to 3 hours	Up to 3 hours
Tablet camera	5MP front, 8MP rear	5MP front, 8MP rear
Dimensions	260 x 215 x 60 mm	260 x 215 x 60 mm
Connectivity	WiFi 6: IEEE 802.11 a/b/g/n/ac/ax Bluetooth 5.0	WiFi 6: IEEE 802.11 a/b/g/n/ac/ax Bluetooth 5.0
Mounting plate	Rehadapt and Daessy	Rehadapt and Daessy
Other features	Screen protector Handle Strap Speaker Tough case Keypad ready USB C to A adapter	Screen protector Handle Strap Speaker Tough case Keypad ready USB C to A adapter