



A slim and lightweight communication aid for touch access, powered by Grid for iPad.

Available in three different sizes, Talk Pad offers users a powerful and portable device with a long battery life and loud speaker.

Tablet	iPad Mini 8.3" 6th generation	iPad 10.2" 8th generation	iPad Pro 12.9" 4th generation
Operating system	iPadOS 15	iPadOS 15	iPadOS 15
Software	Grid for iPad	Grid for iPad	Grid for iPad
Processor	Hexa-core (2× 2.93 GHz Avalanche and 4× 1.82 GHz Blizzard)	Hexa-core (2x2.5 GHz Vortex + 4x1.6 GHz Tempest)	M12 coprocessor Octa-core (4x2.5 GHz Vortex + 4x1.6 GHz Tempest)
Memory / RAM	4GB	3GB	6GB
Hard drive	64GB	32GB	128GB
CPU	A15 Bionic chip with 64-bit architecture	A12 Bionic chip with 64-bit architecture	A12Z Bionic chip with 64-bit architecture
GPU	Apple GPU (5-core graphics)	Apple GPU (4-core graphics)	Apple GPU (8-core graphics)
Display	2266x1488-pixel (326 PPI) Antiglare screen protector	2160x1620-pixel (264 PPI) Antiglare screen protector	2732×2048-pixel (264 PPI) Antiglare screen protector
Access	Touch screen Bluetooth switch	Touch screen Bluetooth switch	Touch screen Bluetooth switch
Weight	815g	1123g	1446g
Mounting	-	Optional Rehadapt and Daessy (92g)	Optional Rehadapt and Daessy (92g)
Battery life	10 hours	10 hours	10 hours
Tablet camera	7MP front, 8MP rear	1.2MP front, 8MP rear	7MP front, 10MP (ultrawide), 12MP rear (wide)
Connectivity	Wi-Fi 6 802.11 ax, Bluetooth 5	Wi-Fi 802.11 a/b/g/n/ac, Bluetooth 4.2	Wi-Fi 802.11 a/b/g/n/ac, Bluetooth 5
Other features	Handle Adjustable stand Shoulder strap Laptop bag Integrated speaker Tough case Keypguard ready Drop tested to 1m		