

# Grid Pad 10s



Grid Pad 10s is a dedicated and multi-access device designed for users who want a portable and rugged communication aid.

The high-quality engineering and unique set of features bring you a versatile device that is perfect for getting your message across quickly and using all the features within Grid.

Operating system	Windows 10 Pro
Software	Grid 3
Processor	Intel® Core™ i5-7Y54, 1.2GHz Base, 3.2GHz Turbo
Memory / RAM	8GB
Hard drive	256GB SSD
Display	10.1" 1920 x 1200 daylight bright display with toughened glass and anti-glare finish, 650 NITS brightness
Second screen	3.9" 480 x 128 second screen with toughened glass and anti-glare finish
Weight	1.45kg (1.5kg with mounting plate) / 3.2lb (3.3lb with mounting plate)
Battery life	8 hours typical use
Charge time	5 hours
Tablet camera	2MP front, 5MP rear
Dimensions	260 x 189 x 60 mm
IP rating	IP54 (dust and splash resistant)
Access	2 x 3.5mm switch ports, 2 x USB 3.0 port, 3.5mm headphone port, touch screen
Environment control	GEWA Infrared
Connectivity	Dual band wireless (AC 8265), Bluetooth 4.2 LE
Mounting plate	Daessy and Rehadapt
Warranty	3 years SmartCare (UK only)
Other features	Antimicrobial case and screen, drop tested to 1m, shoulder strap, flip-out stand and remote power button (with wireless switch access)
Optional extras	Tough cover with handle, keyguard